




# Histopathology



Dr James Pease [j.pease@imperial.ac.uk](mailto:j.pease@imperial.ac.uk)

# Session Plan




## Part 1

- Histopathology
  - Biopsies
  - Resection specimens
  - Frozen sections
- Cytopathology
  - Smears
  - Fine needle aspirates

## Part 2

- Conjugation of antibodies
  - Enzymes
  - Fluorescent markers
  - Magnetic beads
- Use of antibodies in diagnosis
  - Immunohistochemistry
  - ELISA
  - Flow cytometry

# Session Plan



## Part 1

- Histopathology
  - Biopsies
  - Resection specimens
  - Frozen sections
- Cytopathology
  - Smears
  - Fine needle aspirates

## Part 2

- Conjugation of antibodies
  - Enzymes
  - Fluorescent markers
  - Magnetic beads
- Use of antibodies in diagnosis
  - Immunohistochemistry
  - ELISA
  - Flow cytometry




# Histopathology v cytopathology

A **histopathologist** deals with **tissues**.

She/he will examine sections, noting the architecture of the tissue and asking what it tells us about a particular condition.

A **cytopathologist** deals with **cells**.

They are often are the individuals responsible for taking the cells from the patient, preparing them for examination and then delivering their expert diagnosis on the cell sample



# Histopathology

A **histopathologist** deals with **tissues**.

Architecture? Likely medical condition?

Useful for diagnosis and to determine the efficacy of a treatment.

A histopathologist works a variety of tissue samples:

- Biopsies
- Resection specimens
- Frozen sections
- Post-mortems

Real-time examination can inform the surgical process.

# Biopsies



A biopsy embedded in paraffin wax ready to be cut into sections

Biopsies are **small sections** of tissues removed from a patient.

Preserved in a formalin solution which cross-linking proteins and **fixes** the tissue.

Embedded in **paraffin wax** to allow very thin sections (2-3 $\mu$ m thick) to be using a microtome.


Mounted on a microscope slide for further preparation prior to analysis.

Is the tissue normal? Is the tissue inflamed? What is the likely cause? (**1-POM-1-9 Inflammation**)

Is the tissue cancerous and, if it is, what type of cancer is it? (**1-POM-1-7 Cancer**)



# Resection specimens



Resections are taken from tissue that has been removed as part of a surgical procedure and can be processed as for a biopsy.

Resections are used primarily to look at the stage the disease (**1-POM-1-7 Cancer**).

Has the cancer penetrated the bowel wall, spread to the lymph nodes or the liver?

Has all the cancerous tissue been removed or is chemotherapy also needed?




Figure: Resection tissue taken for pathological analysis

# Frozen sections




Figure: Preparation of a frozen section

Frozen sections - taken during surgical procedures and examined in **real time** during the operation.

The freshly taken tissue is frozen by a **cryostat**, cut, mounted on slides and stained as for biopsies.


Provides diagnosis in minutes which can be relayed back to the surgeon to inform the surgery.

For example:

- Is the tissue cancerous?
- Has all the cancerous tissue been removed?
- Is there another pathological process going on?



# Timescales



Key considerations are the relative timescales for a result from the histopathology lab to reach the clinician:

- Frozen section: **30 minutes**
- Biopsies: **2-3 days**
- Resection specimen: **5-7 days**

# Cytopathology




**Human pap smear stained with H&E and showing the presence of chlamydia in the vacuoles.**

Source; National Cancer Institute.

**1-POM-1-10: Histopathology and cytopathology:** Recall situations where histopathology and cytopathology might be used as a diagnostic method, and summarise the main steps involved in processing specimens.

# Fine needle aspirates



A patient undergoing a fine needle aspiration.

## Advantages


Can penetrate relatively inaccessible tissues e.g. thyroid nodule and assess a condition without surgery.

## Disadvantages


Looks at cells in isolation – no information regarding tissue architecture.



# Histopathology and cytopathology are a powerful combination




Skin biopsy taken from an individual with Kaposi's sarcoma.



CD31 immunostaining of a skin biopsy.

# Histopathology and cytopathology are a powerful combination




Diagnosis of reactive lymphadenopathy

Cytopathological analysis of a fine needle aspirate



# Session Plan





## Part 1

- Histopathology
  - Biopsies
  - Resection specimens
  - Frozen sections
- Cytopathology
  - Smears
  - Fine needle aspirates

## Part 2

- Conjugation of antibodies
  - Enzymes
  - Fluorescent markers
  - Magnetic beads
- Use of antibodies in diagnosis
  - Immunohistochemistry
  - ELISA
  - Flow cytometry

# Antibodies –versatile tools for diagnosis





The immunology labs offer 86 different tests of which 53 are used to detect the levels of specific antibodies circulating in patients.

e.g. ds DNA

- systemic lupus erythematosus (SLE)
- Sjögren's syndrome
- rheumatoid arthritis

# Antibody conjugates



## Conjugations


**Enzymes:** e.g. peroxidase, alkaline phosphatase.

**Fluorescent probes.** These can allow the rapid measurement of the levels of molecules within a sample.


**Magnetic beads:** e.g. the purification of cell types.

**Drugs:** e.g. Kadcyra, an anti-**HER2** antibody linked to the cytotoxic chemical **emtansine**.

# Antibodies as diagnostic tools




direct detection




indirect detection

# Use of manufactured antibodies




- Blood group serology (**Haem -Blood transfusion**).
- Immunoassays  
e.g. detection of hormones  
circulating antibodies/antigens.
- Immunodiagnosis  
e.g. infectious diseases  
antibody levels  
IgE (**Immno - Hypersensitivity**).

kidney biopsy



IgG-antigen



Complement  
C3



# ELISA (Enzyme linked immunosorbent assay)




Figure A.6 Janeway's Immunobiology, 8ed. (© Garland Science 2012)

ELISA

Enzyme  
Linked  
ImmunoSorbent  
Assay

- Clinical samples (e.g. adheres to a plastic plate)
- Probe with specific antibody raised against the molecule of interest.
- Enzyme conjugation generates a coloured product.
- Reference to standard curve (**MBC - Haemoglobin**) determine precise concentrations of the molecule in the sample.

# Flow cytometry





Figure A.26 part 1 of 2 Janeway's Immunobiology, 8ed. (© Garland Science 2012)

Cells labelled with **differently** conjugated Abs

Run as a stream of single cells through a laser beam (s).

Colour of light emitted and the forward or side scatter of the laser beam denotes the **identity** of the cell surface molecules expressed and the **size** and **granularity** of the cells.

- anti-CD3<sup>+</sup> T cells – pan T cell marker
- anti-CD4<sup>+</sup> T cells – T helper/cells
- anti-CD8<sup>+</sup> T cells – cytotoxic T cells
- anti-CD19<sup>+</sup> B cells
- anti-CD56<sup>+</sup> Natural Killer (NK) Cells



With many thanks to Dr Keith Gould  
and Professor Robert Goldin